



## Canespro Fungal Nail Treatment Set

### Uses

Removes the infected parts of the nail in 2-3 weeks.

### Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.  
Step 2: Cover the treated nail with the plaster provided and keep covered for 24 hours.

Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.  
Repeat steps 1, 2 and 3 daily for 2-3 weeks.

### Format and Ingredients

Treatment set containing:  
10g tube Canespro Urea 40% ointment, 22 waterproof plasters, 1 plastic nail scraper,

**Category:**  
Footcare

**Manufacturer:**  
BAYER PLC

**Pack size:** 10g

**RRP:** £31.49

**Legal Status:** MD