

# Counter Intelligence PLUS | 2025

The training guide for pharmacy assistants



## Replens Long Lasting Vaginal Moisturiser x 6 Pre-filled Applicators

### Uses

Long lasting vaginal moisturiser. Helps replenish internal vaginal moisture, alleviate vaginal dryness symptoms such as discomfort, itch, irritation and provide more comfortable intimacy. One application provides up to 3 days' relief.

### Dosage

Adults: apply once every 3 days. Remove the applicator from the sealed wrapper and grip firmly by the thick end. Shake down like a thermometer to ensure that the contents are at the thin end. Bend the thin pointed end of the applicator back and forth repeatedly. Snap off the thin pointed end of the applicator and discard the break-away tip. The applicator may be inserted into the vagina while you are in a sitting position or when lying on your back with knees bent. Gently insert the thin end well into the vagina as deeply as it will go comfortably. Squeeze the thick end of the applicator firmly to deposit gel.

### Format and Ingredients

Pre-filled applicator

**Manufacturer:**  
CHURCH & DWIGHT UK LTD

**Pack size:** 6

**RRP:** £10.62

**Legal Status:** MD

A division of  
 GROUPE  
EUROCOM

